Dri-STAT® Enzymatic Bilirubin Reagent

TABLE OF CONTENTS

English	1
Français	3
Deutsch	5
Italiano	8
Español	10
Português	12
Dansk	14
Svenska	17
Norsk	19
Ελληνικά	21
日本語	23
中文	25
Lietuviškai	27
Magyar	30
Polski	32
Türkçe	34
REFERENCES	37

Dri-STAT® Enzymatic Bilirubin Reagent

REF 683802

For *In Vitro* Diagnostic Use

INTENDED USE

Dri-STAT Enzymatic Bilirubin Reagent is intended for use in **in vitro** diagnostic determinations of total bilirubin in serum and plasma.

SUMMARY

The measurement of total bilirubin has been accomplished by numerous methods, most of which are based on the reaction of bilirubin with diazotized sulfanilic acid to produce a blue pigment.^{1,2} Recently, an enzyme has been described which catalyzes the oxidation of bilirubin and produces a loss in the characteristic yellow absorption band^{3,4}. Dri-STAT Enzymatic Bilirubin Reagent employs that enzyme, bilirubin oxidase. In the assay medium provided, total bilirubin can be determined by measuring the decrease in absorbance in the range 405 nm to 465 nm.

CLINICAL SIGNIFICANCE

The measurement of total serum or plasma bilirubin is of value in the diagnosis of hemolytic disorders, biliary obstruction, hepatitis and cirrhosis.

PREPARATION

These instructions apply to reagent preparation for use in the manual procedure presented below. Application data sheets for specific automated instruments may indicate a different procedure which should be followed.

The buffer is provided ready for use. Reconstitution of the Enzyme Reagent should be accomplished by the careful addition of distilled water. The amount of water to be added is specified on the vial label. Mix gently to obtain complete solution. **DO NOT SHAKE TO DISSOLVE.**

STORAGE AND STABILITY

The reagent is stable until the expiration date stated on the label when stored at 2 °C to 8 °C. Reconstituted Enzyme Reagent is stable for eight hours at 20 °C to 25 °C, seven days at 2 °C to 8 °C, or thirty days at -15 °C to -20 °C. The Enzyme Reagent may be frozen only once. Do not thaw at temperatures above 30 °C or allow to thaw without gentle mixing.

SPECIMEN HANDLING

Test specimens can be serum or plasma, collected in a suitable manner. EDTA is the anticoagulant recommended for plasma specimens. Total bilirubin in clear specimens has been reported stable for several hours at room temperature and for up to seven days when refrigerated and protected from light.⁵

EXPECTED VALUES

0.1 to 1.0 mg/dL

It is recommended that each laboratory establish the normal range for the area in which it is located. The values quoted above were derived from 105 donors.

PERFORMANCE

Precision

The within-day precision (± 1 S.D.) of twenty replicates was found to be 3.8%, 0.9%, 0.8%, and 0.9% for mean values of 0.78, 2.23, 3.87, and 7.80 mg/dL, respectively. The day-to-day precision (± 1 S.D.) over a twenty-day period was found to be 4.3%, 0.9%, 1.4%, and 0.9% for mean values of 1.3, 5.7, 10.3, and 18.3 mg/dL, respectively.

Sensitivity

Using the procedure described above, the sensitivity of the assay is such that $\Delta A = 0.028$ is approximately equivalent to 1 mg/dL bilirubin.

Dynamic Range

When used in the manual procedure above, this reagent provides results which are linear to 35 mg/dL. Protocols for automated instruments may call for differences in samples and reagent volumes which may affect the dynamic range. Consult the appropriate application data sheet.

Accuracy

Results obtained with Dri-STAT Enzymatic Bilirubin Reagent in the assay of normal and abnormal sera ($n = 50$) showed a correlation (r) = 0.9997, slope = 1.045, and intercept = 0.18 mg/dL when compared with a widely used Jendrassik method.



LIMITATIONS

In the presence of hemoglobin (100 mg/dL), bilirubin values are reduced by 5%.

A list of drugs and other substances which alter levels of total bilirubin or interfere with its determination has been compiled by Young et al.⁶

PROCEDURE

The reaction sequence employed in the assay of bilirubin is as follows:

The total absorption change in the range 405 nm to 465 nm is proportional to the concentration of bilirubin in the sample.

REAGENTS

Active Ingredients	Concentration in Reagent as Formulated
Enzyme Reagent Bilirubin oxidase (microbial)	25000 U/l

The buffer provided contains 100 mmol/L tris (hydroxymethyl) amino methane, pH 8.2 ± 0.2 .

In addition to Dri-STAT Enzymatic Bilirubin Reagent, performance of the test requires a bilirubin calibrator and devices capable of accurate delivery of samples, standards and reagent, as well as accurate measurement of absorbance.

These instructions are intended for use in tests employing manual instrumentation. Specific instructions for reagent use on selected automated instruments are available on request.

System Parameters

Wavelength	465 nm
Incubation temperature	37 °C
Mode	Absorbance
Absorbance Range	0 to 2 A
Cuvette pathlength	1.0 cm
Reaction time	5 minutes
Sample volume	0.05 mL
Buffer volume	1.00 mL
Enzyme reagent volume	0.08 mL

1. Add 1.0 mL buffer to a clean cuvette. Place the cuvette in the spectrophotometer and adjust absorbance reading to zero.
2. Add 0.05 mL of calibrator or sample and mix well.
3. Record initial absorbances of calibrator and samples (A0)

4. Add 0.08 mL Enzyme Reagent, mix well, and incubate at 37 °C for at least five minutes.
5. Record final absorbances of calibrator and samples (A_f).
6. Calculate $\Delta A = A_0 - A_f$ for calibrator (ΔA_c) and samples (ΔA_s).

CALCULATION

Compute total bilirubin levels as follows:

$$\text{Total Bilirubin (mg/dL)} = \frac{\Delta A_s}{\Delta A_c} \times C^*$$

* C is the value in mg/dL of the calibrator used.

NOTES

1. To avoid drift, final absorbances (A_f) should be read no longer than one hour after Enzyme Reagent addition.
2. If an unusually high total bilirubin level is observed, (> 35 mg/dL), dilute the serum with saline and reassay. Multiply results by dilution factor. A dilution of 3, i.e., 0.1 mL serum + 0.2 mL saline is recommended.
3. This assay may be performed at room temperature by increasing the incubation period to ten minutes.

QUALITY CONTROL PROCEDURES

It is recommended that normal and abnormal controls, or calibrators, be included in each set of assays. Commercially available control material with established data confirmed by a reference laboratory may be used for quality control. The assigned value of the control material must be confirmed by this methodology.

Factors other than reagent which might affect the performance of this test include proper instrument function, cleanliness of glassware and cuvettes, and accuracy of sample or reagent pipetting.

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, California 92821
Printed in U.S.A.